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
SEROTYPING LEGIONELLA PNEUMOPHILA FROM IDEXX QUANTI-TRAY™/ LEGIOLERT®
BY LATEX AGGLUTINATION OR DIRECT FLUORESCENT ANTIBODY MICROSCOPY

1. Scope and Application

1.1. This method is intended for use in serotyping Legionella pneumophila isolates directly from positive wells observed in Quanti-Tray/Legiolert.

1.2. Since not all Legionella pneumophila isolates in Legiolert perform equivalently in serotyping protocols, multiple options are provided.

2. Definitions

2.1. Latex agglutination – a method to detect specific antigens in a sample. Microbeads of latex coated with antigen-specific antibodies bind to antigens in a positive sample and agglutinate, signifying the presence of the antigen.

2.2. Direct fluorescent antibody (DFA) microscopy – a method to detect specific antigens in a sample. Monoclonal antibodies directed against Legionella antigens are conjugated to a fluorochrome, bound to slide-fixed organisms, and imaged using fluorescence microscopy following washing and removal of unbound conjugate.

2.3. Quanti-Tray/Legiolert – a most probable number quantification device used in conjunction with Legiolert growth reagent for the detection of Legionella pneumophila in water samples.

3. Interferences

3.1. No interferences in potable water.

4. Safety

4.1. The analyst/technician must know and observe the normal safety procedures required in a microbiology laboratory preparing, using, and disposing of samples, reagents and materials, and while operating sterilizing equipment.

4.2. Mouth pipetting is prohibited.

5. Equipment and Supplies

5.1. Pipettors capable of delivering 5-1000 µL, sterile pipettor tips, sterile loops.

5.2. Sterile razors or scalpels, disposable alcohol wipes.

5.3. Microfuge tubes with \geq 1.5 mL capacity.

5.4. Microfuge.

6. Reagents

6.1. Latex agglutination reagent kits (Oxoid or comparable); store as per manufacturers guidelines.

6.2. Direct fluorescent antibody kits; store as per manufacturers guidelines.

6.3. Buffered charcoal yeast extract (BCYE) petri plate media; store at 2-10°C.

6.4. Sterile, non-buffered, oxidant-free water.

7. Quality Control

7.1. Quality control of latex agglutination or DFA kits should be conducted as per manufacturers guidelines.

8. Serotyping procedures

8.1. To gain access to a positive Quanti-Tray/Legiolert well, use a scalpel or razor to open the Quanti-Tray from the paper side:

8.1.1. Using a sterile alcohol wipe, clean the paper side of the Quanti-Tray/Legiolert corresponding to the well to be sampled.

8.1.2. Use the same alcohol wipe to sterilize the scalpel or razor.

8.1.3. Cut a slit in the paper side of the tray to open up the well.

8.1.4. Extract the desired volume using a pipettor with disposable tips.

8.1.5. Alternatively, a pipet tip may be used to directly puncture the paper, though the incidence of splattering is greatly increased. Therefore, it is advised that this technique be performed in a biosafety cabinet.

8.2. Procedure Option 1 – Culturing isolates from Quanti-Tray/Legiolert before serotyping.

8.2.1. After opening the target positive well, extract 5 uL of culture and deposit to a BCYE plate.

8.2.2. Streak culture for isolation on BCYE and incubate 2-4 days at 36°C.

8.2.3. Perform serotyping with isolated colonies as per manufacturer's instructions.

8.3. Procedure Option 2 – Serotyping directly out of the Quanti-Tray/Legiolert without additional culturing.

8.3.1. Option 2A – Latex agglutination directly from the well:

8.3.1.1. Add 50 µL of each latex agglutination reagent directly to a unique area of the test surface provided in the kit.

8.3.1.2. Extract 50 µL of material from the opened positive well and add to the latex agglutination reagent.

8.3.1.3. Gently mix up and down in the pipettor 2-4 times, then swirl to observe agglutination as per the manufacturer's instructions.

NOTE: Some isolates detected in Quanti-Tray/Legiolert may have a lower final cell concentration and/or a weak agglutination phenotype. If this is the case, then a user should follow option 2B for concentrating the cells prior to a serotyping test.

8.3.2. Option 2B – Concentrating the cells prior to serotyping:

8.3.2.1. Extract 1.5 mL of culture from the opened positive well and put in a microfuge tube.

NOTE: The goal for this procedure is to produce a 10-fold concentrate, which requires ≥ 1.5 mL of culture. Quanti-Tray/ Legiolert small wells hold ~200 µL of culture and, therefore may be limiting if not pooled, making large wells the primary source of test material. If no large wells are positive, then follow either Option 1 or Option 2a.

8.3.2.2. Spin in a microfuge at ~ 5000 rpm for 1 minute.

8.3.2.3. Remove the cell-free supernatant using a pipettor, taking care not to dislodge the cell pellet.

8.3.2.4. Resuspend the pellet in 150 µL of deionized water.

8.3.2.5. For latex agglutination:

8.3.2.5.1. Add 50 µL of each latex agglutination reagent directly to a unique area of the test

surface provided in the kit.

8.3.2.5.2. Add 50 μL of the concentrate to each agglutination spot and gently mix up and down in the pipettor 2-4 times, then swirl to observe agglutination.

8.3.2.6. For DFA:

8.3.2.6.1. Proceed with manufacturer's instructions using concentrate as the test material.

9. Interpretation and Calculations

9.1. Interpret all serotyping results as per manufacturer's instructions.

10. Waste Management

10.1. It is the laboratory's responsibility to comply with all federal, state and local regulations governing waste management, particularly the biohazard and hazardous identification rules and land disposal regulations. Compliance with all sewage discharge permits and regulations is also required.

10.2. Samples, reference materials and equipment known or suspected to have viable bacteria attached or contained must be sterilized prior to disposal.

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17 Copyright 2016 IDEXX Laboratories