## Coilert\*-18 Test Kit

**Introduction**

Coilert-18 either simultaneously detects total coliforms and *E. coli*; or fecal coliforms in water. It is based on IDEXX's patented Defined Substrate Technology\* (DST\*). When total or fecal coliforms metabolize Coilert-18's nutrient-indicator, ONPG, the sample turns yellow. When *E. coli* metabolize Coilert-18's nutrient-indicator, MUG, the sample also fluoresces. Coilert-18 can simultaneously detect these bacteria at 1 cfu/100 mL within 18 hours even with as many as 2 million heterotrophic bacteria per 100 mL present.

**Storage**

Store at 2-25°C away from light.

**Presence/Absence (P/A) Procedure**

- Add contents of one pack to a 100 mL sample in a sterile, transparent, nonfluorescing vessel.
- Cap vessel and shake.
- If sample is not already at 33–38°C, then place vessel in a 35°C waterbath for 20 minutes
- or, alternatively, a 44.5°C waterbath for 7–10 minutes.
- Incubate at 35±0.5°C for the remainder of the 18 hours.
- Read results according to Result Interpretation table below.

**Quantı-Tray\* Enumeration Procedure**

- Add contents of one pack to a 100 mL room temperature water sample in a sterile vessel.
- Cap vessel and shake until dissolved.
- Pour sample/reagent mixture into a Quanti-Tray or Quanti-Tray/2000 and seal in an IDEXX Quanti-Tray Sealer.
- Place the sealed tray in a 35±0.5°C (or 44.5±0.2°C for fecal coliforms) incubator for 18 hours (prewarming to 35°C is not required). For incubation in a water bath, submerge the Quanti-Tray, as is, below the water level using a weighted ring.
- Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

**Result Interpretation**

<b>Appearance</b>	<b>Result</b>
Less yellow than the comparator <sup>1</sup> when incubated at 35±0.5°C or 44.5±0.2°C	Negative for total coliforms and <i>E. coli</i> ; Negative for fecal coliforms
Yellow equal to or greater than the comparator when incubated at 35±0.5°C	Positive for total coliforms
Yellow equal to or greater than the comparator when incubated at 44.5±0.2°C	Positive for fecal coliforms
Yellow and fluorescence equal to or greater than the comparator when incubated at 35±0.5°C	Positive for <i>E. coli</i>

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Coilert-18 results are definitive at 18–22 hours. In addition, positives for both total coliforms and *E. coli* or fecal coliforms observed before 18 hours and negatives observed after 22 hours are also valid.

**Procedural Notes**

- A slight tinge may be observed when Coilert-18 is added to the sample.
- If excess foam causes problems while using Quanti-Tray, you may choose to use IDEXX Antifoam Solution (Catalog# WAFDB) or IDEXX 120 ml vessels with Antifoam (Catalog# WV120SBAF-200).
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures. For example, samples run in Europe should be incubated at 36±2°C for 18-22 hours.
- When following AFNOR validated method for drinking water testing, place the sealed tray in a 36±2°C incubator for 18 hours (prewarming to 36°C is not required).
- Coilert-18 can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*<sup>2</sup> MPN tables should be used to find Most Probable Numbers (MPNs).
- If a water sample has some background color, compare inoculated Coilert-18 sample to a control blank of the same water sample.
- Coilert-18 can be used for *E. coli* detection (but not coliforms) in marine water. Samples must be diluted at least tenfold. Multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- Coilert-18 is a primary water test. Coilert-18 performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- In samples with excessive chlorine, a blue flash may be seen when adding Coilert-18. If this is seen, consider sample invalid and discontinue testing.
- Aseptic technique should always be followed when using Coilert-18. Dispose of in accordance with Good Laboratory Practices.

**Quality Control Procedures**

The following quality control procedure is recommended for each lot of Coilert-18 when used for total coliform and *E. coli* testing at 35±0.5°C:

- Fill three sterile vessels with 100 mL sterile water and inoculate one with an *E. coli*, one with a total coliform and one with a non-coliform using one of the following products:
  - Quanti-Cult™<sup>3</sup> *E. coli*, *Klebsiella pneumoniae* and *Pseudomonas aeruginosa*
  - A sterile loop of ATCC<sup>4</sup> strains, *E. coli* ATCC 25922 or 11775, *Klebsiella pneumoniae* ATCC 31488 and *Pseudomonas aeruginosa* ATCC 10145 or 27853
- Follow the P/A Procedure or Quanti-Tray Enumeration procedure above.
- Results should match the Result Interpretation table above.

**The following quality control procedure is recommended for each lot of Coilert-18 when used for fecal coliform testing at 44.5±0.2°C:**

- Fill two sterile vessels with 100 mL sterile water and inoculate one with a fecal coliform (*E. coli*) and one with a non-fecal coliform (*Pseudomonas aeruginosa*) using the following strains:
  - A sterile loop of ATCC strains, *E. coli* ATCC 25922 or 11775 and *Pseudomonas aeruginosa* ATCC 10145 or 27853
- Follow the Quanti-Tray Enumeration Procedure above.
- Results should match the Result Interpretation table above.

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<sup>[1]</sup> IDEXX P/A Comparator, catalog # WP104; Quanti-Tray Comparator #WQTC, or Quanti-Tray/2000 Comparator #WQ2KC

<sup>[2]</sup> Eaton AD, Ciesceri LS, Greenberg AE. Standard Methods for the Examination of Water and Wastewater. American Public Health Association, 1998, Washington, DC.

<sup>[3]</sup> Quanti-Cult™ Kulturen - IDEXX catalog # WKIT-1001

<sup>[4]</sup> American Type Culture Collection 1-800-638-6597

<sup>\*</sup>Coilert, Defined Substrate Technology, DST and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. Quanti-Cult is a trademark of Remel Inc.

<sup>©</sup> 2011 IDEXX Laboratories, Inc. All rights reserved.

<sup>Manufactured under one or more of the following U.S. patents:</sup> 5,610,029; 5,518,892; 5,620,895; and 5,753,456. Other U.S. and/or foreign patents issued or pending.

<sup>1</sup>Comparator P/A IDEXX, réf. n° WP104; Comparator Quanti-Tray n° WQTC or Quanti-Tray/2000 Comparator n° WQ2KC

<sup>2</sup>Eaton AD, Ciesceri LS, Greenberg AE. Standard Methods for the Examination of Water and Wastewater. American Public Health Association, 1998, Washington, DC.

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# Colilert<sup>\*</sup>-18

 Colilert<sup>\*</sup>-18/Quantl<sup>\*</sup>-Tray for the quantification of E.Coli and coliforms in drinking water. The AFNOR VALIDATION mark is granted to IDEXX for the following method:

idexx.com/water
afnor-validation.org

Certified by AFNOR Certification
ALTERNATIVE METHOD OF ANALYSIS
IDX 33/01 – 11/09
Westbrook, Maine 04092 USA
One IDEXX Drive
IDEXX Laboratories, Inc.

Australia: 1 800 655 978
Japan: + 81-422-71-5921
China: + 86-21-6161-279328
UK: + 44 (0) 1638 676800
Europe: 00800 4339 9111
North/South America: 207-556-4496/1-800-321-0207

**For IDEXX Technical Support, please call:**

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© 2011 IDEXX Laboratories, Inc. Todos los derechos reservados.